


# LED Upgrade

It's light years ahead!



**a dec**<sup>®</sup>



# Upgrade to Brilliance

**Bright, white, and true.** Replacing your current A-dec light head with the new A-dec LED is quick, easy, and entirely hassle-free. Upgrade to 25% more illuminance while consuming one-fifth of the power needed for halogen bulbs.

**Increase comfort, performance, and reduce your cost of ownership.**

- Reduces shadowing, improves contrast and eliminates eye strain
- See colors accurately for better tissue diagnosis
- Gain stunning white light with a color rendering index of 90 at 5,000K
- Optimize ergonomics for the entire dental team
- Position with three axes of rotation and 540° of horizontal motion
- Illuminate without curing with a cure-safe mode at 25,000 lux
- Lower energy costs by 80% and never replace a bulb again
- Treat patients in comfort—without radiating heat or noise



**Am I eligible for an upgrade?** As long as you have a working flexarm on an A-dec light manufactured after 1993, you can upgrade it to an A-dec LED.\*

**Will my new LED be compatible with my A-dec touchpad?** Yes. A-dec touchpads utilizing DCS (Data Communication System) give you the choice to operate the light using the touchpad or local switching. With A-dec DCS, you can power on/off, select from three intensities, activate cure-safe mode, and support auto on/off with chair presets. Without DCS, your new A-dec LED light will have similar control via local switching only.

**You can Upgrade to Brilliance now for RRP (inc. GST) \$3,115**

**How can I learn more?** Visit [a-dec.com/LED](http://a-dec.com/LED). Contact your Authorized A-dec Dealer or phone A-dec on 1800 225 010.

## Ordering Information

Description	Part Number	RRP Price (inc. GST)
A-dec LED light head upgrade in White	28.1924.00	\$3,115
A-dec LED light head upgrade in Surf	28.1925.00	\$3,115

\* Compatibility exceptions apply for the A-dec Performer Light and some A-dec 300 lights. For details, see your Authorized A-dec Dealer. Installation costs may apply.



Chairs  
Delivery Systems  
Lights  
Monitor Mounts  
Cabinets  
Handpieces  
Maintenance  
Sterilisation  
Imaging

**For more information**  
Phone: 1800 225 010

Email: [a-dec@a-dec.com.au](mailto:a-dec@a-dec.com.au)  
Visit: [www.a-dec.com.au](http://www.a-dec.com.au)